

Using on-line UV-VIS spectroscopy for rinse monitoring with the new developed flowcell R.V.P. (**R**epeatable **V**ariable **P**ath length)

Cleaning Validation is a critical aspect in pharmaceutical batch manufacturing. Production run campaigns are followed by downtime periods of 'Change-Over' and this can have a considerable impact on production capacity and profitability.



Fig. 1: CV online installation in a pharmaceutical plant.

Normally samples from the wash cycles are measured off-line by HPLC. This processes are time consuming and costly.

Using Process Analytical Technology (PAT) philosophies, is an alternative the installation of an on-line analyser. High sensitivity and fast response time allow the monitoring of cleaning processes in real time. A quick feed back is guaranteed that allows optimum cycle times, reducing dead times and saving costs.

On behalf of affiliated global pharmaceutical groups J&M Analytik AG designed around it's new R.V.P.- flow cell both, close fitted or mobile analysers.





Fig. 2 & Fig. 3: R.V.P. flow cell installed on a mobile CIP- analyzer for bypass operation.

In primary or secondary manufacturing plants, production equipment with cleaning in place (CIP) systems has to be thoroughly cleaned by a validated process before the next production cycle can start. An on- line CIP analyzer gives you quick feed back and ensures that your cycle times are optimal.

Applications

Process monitoring and controlling in chemical-, pharmaceutical-, food industries, biotechnology research and production, environmental monitoring, nearly unlimited (...)

Your approach to implement PAT for cleaning in place

| PAT for cleaning process | | QC Laboratory |
|--------------------------|---|---|
| Stage 1 | Use process analyzer(s) to visualize and gain insight in the cleaning process. Optimize the cleaning process. IMS, DART, UV, NIR, ... | Development of innovative analytical technologies to reduce time for method development, validation and sample analysis e.g. UPLC, Challenge process analyzer(s), GMP-analysis of cleaning samples for release of equipment |
| Stage 2 | Routine monitoring of optimized cleaning process. Evaluate performance of PAT tool(s). | GMP-analysis of cleaning samples for release of equipment. |
| Stage 3 | In-line cleaning analysis and immediate release of equipment. | SWOP TEST |
| Stage 4 | Real-time control of cleaning process. | SWOP TEST |



Fig. 3:
Typical mobil CIP process analyzer J&M detector inside



Fig. 4:
Typical laboratory setup (HPLC) J&M detector inside




Fig. 5
Graph shows the change in intensity over time in the entire wavelength range from 190 to 400 nm. The spectra were acquired every 5 sec.




Fig. 6
Graph shows the change in intensity over time for the wavelength at 265 nm.

This example proves the advantages of using J&M instrumentation and software in online analysis applications like reaction monitoring, cleaning validation and blending or mixing procedures.

Please note:
Due to technological progress all specifications can be changed without a further note.